

EYESPRAY

The new type of eyewash

- The eyewash is activated automatically by pulling the cover and letting it fall down.
- The shower is loaded with two spray bottles- each bottle contains a sterile 0,9 % buffered sodium chloride solution.
- The flushing time exceeds 15 min.
- When the shower has been activated both your hands are free to keep your eyelids open while your eyes and face is flushed.
- **Requires no adjustment to the water mains, which eliminates all eventual problems with Legionella bacteria.**
- During transportation of the injured person to hospital separate spray bottles should be used.
- Suitable applications are industries, laboratories, hospitales and all enviroments with limited access to the water mains.
- Each bottle has an expiration date after 36 month from the production date.



If you have any questions please give us a call or send us an e-mail.

